


GCaMP6f

2D Acquisition/Analysis Workflow




Time Series Acquisition

The camera is configured to turn on the light source, take a high resolution image, save it as a TIF (with the timestamp as its filename), turn off the light source, and wait 10 seconds before repeating.

Typically a 10m baseline period is used, 3m of 200 μ M TGOT, then 10m of washout. While increasing the frame rate would increase the number of data points, it also increases accumulated light exposure (photobleaching) and data storage requirements.

Information about the experiment (baseline, drug times, etc.) is stored in `experiment.txt`

Early experiments analyzed ~1hr of data consisting of high resolution TIFs acquired every 500ms. This is approximately 8GB of data. These frames make excellent video, but have far more data points than are required to draw an experimental conclusion.




Manual Selection of ROIs

Although dozens of fluorescent cells are visible, not all visible cells respond to the drug of interest.

ROIs are outlined by hand with FIJI/ImageJ.

- The first ROI is of a region of tissue which contains no cells. It is used as a baseline to subtract-out photobleaching in subsequent ROIs.
- The ROI file is always saved as `RoiSet.zip`
- The *multi-measure* command is used to analyze every ROI of every frame (reported as average pixel intensity) and saved as `Results.xls`

Not every fluorescent cell is a healthy cell. Dead and dying cells become fluorescent due to a loss of calcium regulation capability. Typically the best responding cells are dark initially and become fluorescent upon activation. Healthy quiescent cells will exhibit low fluorescence.




Automated Analysis

If these steps were followed as written, a single experiment will result in a single folder containing 100s of TIFs, `experiment.txt`, `RoiSet.zip`, and `Results.xls`. Analysis is automatic and occurs when the user goes to the web page for this folder.

Invisible steps that occur when you load the web page:

- An R script creates `results_B.xls` which contains ROI data in baseline-subtracted $\Delta F/F$ units.
- An R script creates graphs of ROI data in $\Delta F/F$ units.
- A Python script combines these results with the series of TIFs to create a HTML5-compatible video (`render.mp4`)



GCaMP6f

Fluorescence Quantification Method

Raw Intensity Value
When a ROI is converted to a single value per frame, we refer to it as the raw value (f).
When these values are strung together, it graphs of a single ROI with respect to time.

